

Sterilization Without Compromise

- **Expertise** – Specializing in Sterilization
- **Experience** - Over 90 Years
- **Choice** - Largest Line to Choose from

ALL YOU NEED TO DO IS CHOOSE



Fully Automatic Autoclaves



Elara11 - Class B
Pre and Post Vacuum Autoclaves



Large Capacity Autoclaves



Manual Autoclaves



Compact Autoclaves



Water Distillers



Ultrasonic Cleaners
and Tablets



Autoclave Cleaner

| Model | Chamber Dimensions Ø x L in(mm) | Volume | Overall Dimensions DxWxH in(mm) | Standard Cassette Capacity | Tray Dimensions DxWxH in(mm) | No. of Trays | Standard Unwrapped Cycle Time | Voltage(V) Freq.(Hz) | Power(W) Current(A) | Shipping Weight lbs(kg) |
|-------|---------------------------------------|--------|---------------------------------------|----------------------------|------------------------------------|--------------|-------------------------------|-------------------------|------------------------|----------------------------|
|-------|---------------------------------------|--------|---------------------------------------|----------------------------|------------------------------------|--------------|-------------------------------|-------------------------|------------------------|----------------------------|

Fully Automatic Autoclaves

| | | | | | | | | | | |
|-----------------|--------------------------|-----------------|--|--|-------------------------------------|---|-------------------------------|---------------|--------------|-------------|
| EZ9Plus | 9" x 19.8" (230x504) | 5.2gal 19.8L | 24.8" x 19.1" x 15.1" ^{**} (630x485x384) | 2 Full & 2 Half ^{**} (loaded horizontally) | 16.3" x 6.7" x 0.6" (413x169x16) | 3 | Cold: 30 mins Hot: 24 mins | 120V 50/60 | 1400W 12A | 104 (47) |
| EZ11Plus | 11" x 19.8" (280x504) | 7.5gal 28.5L | 24.8" x 20.9" x 17.3" ^{**} (630x530x440) | 4 Full & 4 Half ^{**} (loaded vertically) | 16.3" x 6.7" x 0.6" (413x169x16) | 5 | Cold: 37 mins Hot: 28 mins | 120V 50/60 | 1400W 12A | 126 (57) |

*Models EZ9Plus and EZ11Plus require a minimum tabletop depth of 22".

| | | | | | | | | | | |
|------------------------------|--------------------------|-----------|--------------------------------------|--|--|---|-------------------------------|----------------|----------------|-------------------------------------|
| EZ9 | 9" x 18" (230 x 470) | 5gal/19L | 21.5" x 20" x 14.4" (545x510x365) | 2 full ^{**} and 2 half ^{**} | 16.3" x 6.7" x 0.8" (415x170x20) | 3 | Cold: 27 mins Hot: 13 mins | 120V 50/60 | 1400W 12A | 95 (43) |
| EZ10 | 10" x 19" (254 x 475) | 6gal/23L | 21.5" x 20" x 14.4" (545x510x365) | 3 full ^{**} and 3 half ^{**} | 16.3" x 6.7" x 0.8" (415x170x20) | 4 | Cold: 30 mins Hot: 14 mins | 120V 50/60 | 1400W 12A | 95 (43) |
| EZ10k | 10" x 19" (254 x 475) | 6gal/23L | 21.5" x 20" x 14.4" (545x510x365) | 3 full ^{**} and 3 half ^{**} | 16.3" x 6.7" x 0.8" (415x170x20) | 4 | Cold: 21 mins Hot: 11 mins | 230V* 50/60 | 2200W 10A | 95 (43) |
| 3870EA with stand | 15" x 30" (384 x 760) | 22gal/85L | 34.5" x 26" x 20.6" (875x665x525) | N/A | 26.6" x 11.3" x 1" (675x286x25) 26.6" x 13.8" x 1" (675x350x25) | 2 | Cold: 32 mins Hot: 23 mins | 230V* 50/60 | 3000W 13.6A | 278 (126) (includes stand) |

Manual Autoclaves

| | | | | | | | | | | |
|---------------------------------|--------------------------|-----------|--------------------------------------|--|--|---|-------------------------------|----------------|----------------|-------------------------------------|
| Valueklave 1730 | 7" x 13" (170 x 340) | 2gal/7.5L | 17.9" x 17.4" x 12" (455x440x305) | 2 half ^{**} | 11.6" x 4.7" x 0.8" (295x120x20) | 3 | Cold: 16 mins Hot: 11 mins | 120V 50/60 | 1350W 11A | 55 (25) |
| 2340M | 9" x 18" (230 x 470) | 5gal/19L | 21.5" x 20" x 14.4" (545x510x365) | 2 full ^{**} and 2 half ^{**} | 16.3" x 6.7" x 0.8" (415x170x20) | 3 | Cold: 27 mins Hot: 13 mins | 120V 50/60 | 1400W 12A | 85 (39) |
| 2540M | 10" x 19" (254 x 475) | 6gal/23L | 21.5" x 20" x 14.4" (545x510x365) | 3 full ^{**} and 3 half ^{**} | 16.3" x 6.7" x 0.8" (415x170x20) | 4 | Cold: 30 mins Hot: 14 mins | 120V 50/60 | 1400W 12A | 85 (39) |
| 2540MK | 10" x 19" (254 x 475) | 6gal/23L | 21.5" x 20" x 14.4" (545x510x365) | 3 full ^{**} and 3 half ^{**} | 16.3" x 6.7" x 0.8" (415x170x20) | 4 | Cold: 21 mins Hot: 11 mins | 230V* 50/60 | 2200W 10A | 85 (39) |
| 3870M with stand | 15" x 30" (384 x 760) | 22gal/85L | 34.5" x 26" x 20.6" (875x665x525) | N/A | 26.6" x 11.3" x 1" (675x286x25) 26.6" x 13.8" x 1" (675x350x25) | 2 | Cold: 32 mins Hot: 23 mins | 230V* 50/60 | 3000W 13.6A | 278 (126) (includes stand) |

Pre and Post Vacuum (Class B) Autoclave

| | | | | | | | | | | |
|----------------|--------------------------|-------------------|--|--|-------------------------------------|---|-------------------------------|----------------|--------------|----------|
| Elara11 | 11" x 19.8" (280x504) | 7.5 gal. 28.5L | 26" x 20.9" x 17.3" ^{**} (660x530x440) | 4 Full & 4 Half ^{**} (loaded vertically) | 16.3" x 6.7" x 0.6" (442x173x19) | 5 | Cold: 25 mins Hot: 18 mins | 230V* 50/60 | 2300W 10A | 171 (78) |
|----------------|--------------------------|-------------------|--|--|-------------------------------------|---|-------------------------------|----------------|--------------|----------|

*Model Elara11 requires a minimum tabletop depth of 24"

**Requires a voltage between 220V and 235V for proper operation. When applicable a buck/boost transformer (0.5KVA) may be required.

** Standard full size cassette is 8"x11" and standard half size cassette is 8"x5.5"



Tuttnauer USA Co., Ltd.

25 Power Drive, Hauppauge, NY 11788 Tel: +800 624 5836, +631 737 4850 Fax: +631 737 0720

E-mail: info@tuttnauerUSA.com, www.tuttnauerUSA.com

Tuttnauer
Your Sterilization & Infection Control Partners